

SPECTROGON

Optical filters • Coatings • Gratings

Date: 2021-11-05

Operator: EBO

Article:

Art no: 711M

LP-0565 nm Ø25.4x1.5 mm

Specifications:

Cuton 5 %: 565 ±20 nm


Tavg: >75 %

Blocking: 200 nm

Measurement conditions:

Temperature: 23 °C

AOI: 0 deg



Measured values:

Cuton 5 % = 574,66 nm...OK

HP1 = 598,84 nm

Tavg = 83,08 %...OK

Slope 1 = 5,27 %

Spectrogon:

M nr: 64497

Calibration traceable to NIST.

Spectrogon products conform in all aspects with the RoHS Directive.