




SPECTROGON

Optical filters • Coatings • Gratings

Article:

NB-1765-035 nm

Specifications:

CWL: +10/-10 nm

HW: +5/-5 nm

Tmin: 70 %

Blocking: 3500 nm

Measurement conditions:

Temperature: 23 °C

AOI: 0 deg

Actual values:

CWL = 1768,6 nm

HW = 35 nm

Cuton 5% = 1731,8 nm

Cutoff 5% = 1805,9 nm

HP1 = 1751,1 nm

HP2 = 1786,1 nm

Tpeak = 78,16 %

Tavg = 73,18 %

Slope 1 = 1,4 %

Slope 2 = 1,5 %

Data collection range: 1000 - 2050 nm

Diagram range: 1500 - 2000 nm

Ref: 4-4297

Calibration traceable to NIST.

Spectrogon products conform in all aspects with the RoHS Directive.